

Anti-Mrps18C Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q9Y3D5

Item Number: CM0021920



Introduction

Rabbit polyclonal antibody against MRPS18C, a mitochondrial ribosomal protein of the small 28S subunit. Validated for WB, IHC-P, IF/ICC and ELISA; reacts with human, mouse, and rat. For mitochondrial translation research.

[Learn More](#)

Parameter	Value
Product Name	MRPS18C Rabbit pAb
Remarks/Alias	CGI-134; S18mt-c; MRPS18-1; mrps18-c; MRP-S18-1; MRP-S18-c; MRPS18C
Species	Human
GeneID (Human)	51023
GeneID	51023
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-142 of human MRPS18C (NP_057151.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y3D5
Protein Weight	16kDa
Shipping	Ice bag