

Anti-Mrps2 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q9Y399

Item Number: CM0015344



Introduction

This rabbit polyclonal antibody recognizes human MRPS2 (S2mt), a component of the mitochondrial ribosome essential for protein synthesis. Validated for WB, IHC-P, and ELISA; also detects mouse MRPS2. Supports mitochondrial translation and disease studies.

[Learn More](#)

Parameter	Value
Product Name	MRPS2 Rabbit pAb
Remarks/Alias	S2mt; CGI-91; MRP-S2; COXPD36; MRPS2
Species	Human
GeneID (human)	51116
GeneID	51116
Immunogen	Recombinant fusion protein containing amino acids 1-296 of human MRPS2 (NP_057118.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Mouse
SWISS	Q9Y399
Protein Weight	33kDa
Shipping	Ice bag