


Anti-Mrps21 Rabbit Polyclonal Antibody For Ihc, Elisa - P82921

Item Number: CM0019196



Introduction

Rabbit polyclonal antibody targeting mitochondrial ribosomal protein MRPS21 (bS21m). Validated for IHC-P and ELISA with reactivity in human, mouse, and rat samples. Supports mitochondrial biology and protein synthesis studies.

[Learn More](#)

Parameter	Value
Product Name	MRPS21 Rabbit pAb
Remarks/Alias	bS21m; MDS016; RPMS21; MRP-S21
Species	Human
Gene ID (Human)	54460
Gene ID	54460
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-87 of human MRPS21(NP_114107.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P82921
Protein Weight	10kDa
Shipping	Ice bag