


Anti-Mrps22 Rabbit Polyclonal Antibody For Wb, Elisa - P82650

Item Number: CM0028658



Introduction

Rabbit polyclonal antibody targeting MRPS22 (mitochondrial ribosomal protein S22) validated for WB and ELISA. Shows cross-reactivity with human, mouse, and rat. ~41 kDa. Suitable for mitochondrial translation and primary mitochondrial disease research.

[Learn More](#)

Parameter	Value
Product Name	MRPS22 Rabbit pAb
Remarks/Alias	GIBT; ODG7; GK002; C3orf5; COXPD5; RPMS22; MRP-S22; MRPS22
Species	Human
GeneID	56945
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-360 of human MRPS22 (NP_064576.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P82650
Protein Weight	41kDa
Shipping	Ice bag