

Anti-Mrps23 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q9Y3D9

Item Number: CM0023807

Introduction

High-quality rabbit polyclonal antibody against human MRPS23, a mitochondrial ribosomal protein implicated in primary mitochondrial disease and COXPD46. Validated for Western blot, IF/ICC, and ELISA applications.

[Learn More](#)

Parameter	Value
Product Name	MRPS23 Rabbit pAb
Remarks/Alias	CGI-138; COXPD46; HSPC329; MRP-S23; MRPS23
Species	Human
GeneID (Human)	51649
GeneID	51649
Immunogen	Recombinant protein corresponding to a sequence within amino acids 1-190 of human MRPS23 (NP_057154.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	Q9Y3D9
Protein Weight	22kDa
Shipping	Ice bag