

Anti-Mrps26 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q9Byn8

Item Number: CM0030745



Introduction

High-quality rabbit polyclonal antibody against human MRPS26, validated in WB, IHC-P, ELISA. Reacts with human, mouse, rat. Ideal for mitochondrial ribosome research and protein expression studies.

[Learn More](#)

Parameter	Value
Product Name	MRPS26 Rabbit pAb
Remarks/Alias	GI008; MRPS13; RPMS13; MRP-S13; MRP-S26; NY-BR-87; C20orf193; dj534B8.3; MRPS26
Species	Human
GeneID (Human)	64949
GeneID	64949
Immunogen	Recombinant protein: Recombinant fusion protein containing a sequence corresponding to amino acids 1-205 of human MRPS26 (NP_110438.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9BYN8
Protein Weight	24kDa
Shipping	Ice bag