

Anti-Mrps28 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q9Y2Q9

Item Number: CM0028363



Introduction

Rabbit polyclonal antibody against human MRPS28 (bS1m), validated for WB, IF/ICC, and ELISA. Cross-reacts with human and mouse. Suitable for mitochondrial translation and disease research.

[Learn More](#)

Parameter	Value
Product Name	MRPS28 Rabbit pAb
Remarks/Alias	MRPS35; COXPD47; HSPC007; MRP-S28; MRP-S35; MRPS28
Species	Human
Gene ID (Human)	28957
Gene ID	28957
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-187 of human MRPS28 (NP_054737.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9Y2Q9
Protein Weight	21kDa
Shipping	Ice bag