


Anti-Mrps30 Polyclonal Antibody For Wb, If/Icc, Ip, Elisa - Q9Np92

Item Number: CM0024892



Introduction

High-quality rabbit polyclonal antibody against human MRPS30 (also known as PAP, PDCD9, S30mt). Validated for WB, IF/ICC, IP, ELISA applications. Reacts with human, mouse, rat mitochondrial ribosomal protein mL65.

[Learn More](#)

Parameter	Value
Product Name	MRPS30 Rabbit pAb
Remarks/Alias	PAP; PDCD9; S30mt; MRP-S30; MRPS30
Species	Human
GeneID (Human)	10884
GeneID	10884
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 190-439 of human MRPS30 (NP_057724.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9NP92
Protein Weight	50kDa
Shipping	Ice bag