


Anti-Mrps31 Rabbit Monoclonal Antibody For Wb, Elisa - Q92665

Item Number: CM0002671



Introduction

High-quality rabbit monoclonal antibody against human MRPS31 (mitochondrial ribosomal protein S31). Validated for WB and ELISA, with cross-reactivity to mouse and rat. Ideal for mitochondrial translation studies.

[Learn More](#)

Parameter	Value
Product Name	MRPS31 Rabbit mAb
Remarks/Alias	S31mt; IMOGN38; MRP-S31; MRPS31
Species	Human
GeneID (Human)	10240
GeneID	10240
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 296-395 of human MRPS31 (Q92665).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q92665
Protein Weight	45kDa
Shipping	Ice bag