

Anti-Mrps31 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q92665

Item Number: CM0016539



Introduction

Rabbit polyclonal antibody targeting mitochondrial ribosomal protein MRPS31 (S31mt, MRP-S31). Validated for Western blot, immunofluorescence, and ELISA in human and rat samples. Detects a band at ~45 kDa.

[Learn More](#)

Parameter	Value
Product Name	MRPS31 Rabbit pAb
Remarks/Alias	S31mt; IMOGN38; MRP-S31; MRPS31
Species	Human
Gene ID (Human)	10240
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 66-395 of human MRPS31 (NP_005821.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Rat
SWISS	Q92665
Protein Weight	45kDa
Shipping	Ice bag