


Anti-Mrps35 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P82673

Item Number: CM0017624



Introduction

High-quality rabbit polyclonal antibody specific to human MRPS35, a mitochondrial ribosomal protein. Validated for WB, IHC-P, IF/ICC, and ELISA. Cross-reacts with mouse. Supports mitochondrial ribosome research.

[Learn More](#)

Parameter	Value
Product Name	MRPS35 Rabbit pAb
Remarks/Alias	MDS023; MRPS28; MRP-S28; HDCMD11P; MRPS35
Species	Human
GeneID (Human)	60488
GeneID	60488
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-323 of human MRPS35 (NP_068593.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	P82673
Protein Weight	37kDa
Shipping	Ice bag