


Anti-Mrps6 Rabbit Polyclonal Antibody For Wb And Elisa - P82932

Item Number: CM0028953



Introduction

High-quality rabbit polyclonal antibody targeting human MRPS6, validated for WB and ELISA. Cross-reacts with mouse and rat. Suitable for studying mitochondrial ribosomal protein S6 in translation and ribosomal biology. 14 kDa protein, P82932.

[Learn More](#)

Parameter	Value
Product Name	MRPS6 Rabbit pAb
Remarks/Alias	S6mt; RPMS6; MRP-S6; C21orf101; MRPS6
Species	Human
GeneID (Human)	64968
GeneID	64968
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-125 of human MRPS6 (NP_115865.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P82932
Protein Weight	14kDa
Shipping	Ice bag