


Anti-Mrps7 Polyclonal Antibody For Wb, Elisa - Q9Y2R9

Item Number: CM0018916



Introduction

KD validated rabbit polyclonal antibody against mitochondrial ribosomal protein S7 (MRPS7), suitable for Western blot and ELISA applications. Targets human MRPS7, protein weight ~28 kDa. Ideal for mitochondrial research and disease studies.

[Learn More](#)

Parameter	Specification
Product Name	[KD Validated] MRPS7 Rabbit pAb
Remarks/Alias	MRPS7; MRP-S; MRP-S7; RP-S7; RPMS7; S7mt; bMRP27a; mitochondrial ribosomal protein S7
Species	Human
GeneID (Human)	51081
GeneID	51081
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 40-140 of human MRPS7 (NP_057055.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q9Y2R9
Protein Weight	28kDa
Shipping	Ice bag