

## Anti-Mrps7 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q9Y2R9

Item Number: CM0023361



### Introduction

Validated for WB, IF/ICC, ELISA; reacts with human, mouse, rat MRPS7 (28kDa, Q9Y2R9). Suitable for mitochondrial ribosome studies and primary mitochondrial disease research.

[Learn More](#)

| Parameter        | Value                                                                                                              |
|------------------|--------------------------------------------------------------------------------------------------------------------|
| Product Name     | MRPS7 Rabbit pAb                                                                                                   |
| Remarks/Alias    | S7mt; MRP-S; RP-S7; RPMS7; MRP-S7; COXPD34; bMRP27a; MRPS7                                                         |
| Species          | Human                                                                                                              |
| GeneID (human)   | 51081                                                                                                              |
| GeneID           | 51081                                                                                                              |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 110-190 of human MRPS7 (NP_057055.2) |
| Source           | Rabbit                                                                                                             |
| Category         | Polyclonal Antibodies                                                                                              |
| Application      | WB, IF/ICC, ELISA                                                                                                  |
| Cross Reactivity | Human, Mouse, Rat                                                                                                  |
| SWISS            | Q9Y2R9                                                                                                             |
| Protein Weight   | 28kDa                                                                                                              |
| Shipping         | Ice bag                                                                                                            |