


Anti-Mrps9 Polyclonal Antibody For Wb, Ihc-P, Elisa - P82933

Item Number: CM0017367



Introduction

Rabbit polyclonal antibody against human MRPS9 (mitochondrial ribosomal protein S9) for WB, IHC-P, ELISA. Cross-reacts with human, mouse, rat. Ideal for mitochondrial translation and ribosome research.

[Learn More](#)

Parameter	Value
Product Name	MRPS9 Rabbit pAb
Remarks/Alias	S9mt; RPMS9; MRP-S9; MRPS9
Species	Human
GeneID (human)	64965
GeneID	64965
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 237-396 of human MRPS9 (NP_872578.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P82933
Protein Weight	46kDa
Shipping	Ice bag