


Anti-Mrrf Monoclonal Antibody For Wb, Elisa - Q96E11

Item Number: CM0007551



Introduction

MRRF Rabbit monoclonal antibody validated for Western blot and ELISA. Detects mitochondrial ribosome-recycling factor in human, mouse, and rat samples. Ideal for mitochondrial translation studies.

[Learn More](#)

Parameter	Value
Product Name	MRRF Rabbit mAb
Remarks/Alias	RRF; MRFF; MTRRF; MRRF
Species	Human
GeneID (Human)	92399
GeneID	92399
Immunogen	Synthetic Peptide
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q96E11
Protein Weight	29kDa
Shipping	Ice bag