

Anti-MS4A14 Rabbit Monoclonal Antibody For Wb, Elisa - Q96Ja4

Item Number: CM0002621

Introduction

Rabbit monoclonal antibody targeting human MS4A14 (Q96JA4), validated for WB and ELISA with rat cross-reactivity. Designed for research on membrane receptor complexes and signal transduction pathways.

[Learn More](#)

Parameter	Value
Product Name	MS4A14 Rabbit mAb
Remarks/Alias	MS4A16; NYD-SP21; MS4A14
Species	Human
GeneID (Human)	84689
GeneID	84689
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 150-280 of human MS4A14 (Q96JA4).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Rat
SWISS	Q96JA4
Protein Weight	77kDa
Shipping	Ice bag