

Anti-Ms4A4A Rabbit Polyclonal Antibody - Ms4A4A

Item Number: CM0022257

Introduction

Rabbit polyclonal antibody targeting mouse Ms4a4a, validated for Western blot and ELISA. High specificity for research on Ms4a4a in mouse samples. Ideal for immunodetection studies.

[Learn More](#)

Parameter	Description
Product Name	Ms4a4a Rabbit pAb
Remarks/Alias	EG666907; Ms4a4a
Host Species	Rabbit
Gene ID (Mouse)	666907
Human Ortholog Gene ID	51338
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 150-236 of mouse Ms4a4a (NP_001297260.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
Shipping	Ice bag