

Anti-MS4A7 Polyclonal Antibody For Wb, Elisa - Q9Gzw8

Item Number: CM0015520

Introduction

Rabbit polyclonal antibody targeting human MS4A7, validated for Western blot and ELISA. Cross-reacts with mouse. Ideal for studies on signal transduction, monocyte differentiation, and tissue expression profiling.

[Learn More](#)

Parameter	Value
Product Name	MS4A7 Rabbit pAb
Remarks/Alias	CFFM4; MS4A8; 4SPAN2; CD20L4; MS4A7
Species	Human
GeneID (Human)	58475
GeneID	58475
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 100-200 of human MS4A7 (NP_067024.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	Q9GZW8
Protein Weight	26kDa
Shipping	Ice bag