

## Anti-MS4A8B Polyclonal Antibody For Wb, Elisa - Q9By19

Item Number: CM0017753

### Introduction

CamelBio offers this anti-MS4A8B polyclonal antibody, validated for WB/ELISA. It targets human MS4A8B, a four-span transmembrane protein expressed in hematopoietic tissues, supporting signal transduction studies.

[Learn More](#)

| Parameter        | Value                                                                                                     |
|------------------|-----------------------------------------------------------------------------------------------------------|
| Product Name     | MS4A8B Rabbit pAb                                                                                         |
| Remarks/Alias    | MS4A4; 4SPAN4; CD20L5; MS4A8B                                                                             |
| Species          | Human                                                                                                     |
| Gene ID (Human)  | 83661                                                                                                     |
| Gene ID          | 83661                                                                                                     |
| Immunogen        | A synthetic peptide corresponding to a sequence within amino acids 100-200 of human MS4A8B (NP_113645.1). |
| Source           | Rabbit                                                                                                    |
| Category         | Polyclonal Antibodies                                                                                     |
| Application      | WB, ELISA                                                                                                 |
| Cross Reactivity | Human                                                                                                     |
| SWISS            | Q9BY19                                                                                                    |
| Protein Weight   | 26kDa                                                                                                     |
| Shipping         | Ice bag                                                                                                   |