


Anti-Msh2 Monoclonal Antibody For Wb, Elisa - P43246

Item Number: CM0008440



Introduction

Rabbit monoclonal antibody targeting human MSH2 (DNA mismatch repair protein Msh2) for western blot (WB) and ELISA applications. Detects human, mouse, and rat MSH2 protein at 105 kDa.

[Learn More](#)

Parameter	Value
Product Name	MSH2 Rabbit mAb
Remarks/Alias	FCC1; COCA1; HNPCC; LCF52; MSH-2; hMSH2; HNPCC1; LYNCH1; MMRC52; MSH2
Species	Human
GeneID (Human)	4436
GeneID	4436
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human MSH2 (P43246).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P43246
Protein Weight	105kDa
Shipping	Ice bag