


Anti-Msh2 Polyclonal Antibody For Wb, If/Icc, Elisa - P43246

Item Number: CM0012954



Introduction

Rabbit polyclonal antibody against human MSH2, validated for WB, IF/ICC, ELISA with cross-reactivity to mouse and rat. Ideal for DNA mismatch repair, cancer research, and cell cycle studies.

[Learn More](#)

Parameter	Value
Product Name	MSH2 Rabbit pAb
Remarks/Alias	FCC1; COCA1; HNPCC; LCFS2; MSH-2; hMSH2; HNPCC1; LYNCH1; MMRC52; MSH2
Species	Human
Gene ID (Human)	4436
Gene ID	4436
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of human MSH2 (NP_000242.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P43246
Protein Weight	105kDa
Shipping	Ice bag