


## Anti-Msh2 Rabbit Polyclonal Antibody (Ko Validated) For Wb, If/Icc, Ip, Elisa - P43246

Item Number: CM0022618



### Introduction

Rabbit polyclonal MSH2 antibody, KO validated, suitable for Western blot, immunofluorescence, immunoprecipitation, and ELISA. Detects human and rat MSH2, a key DNA mismatch repair protein, supporting cancer and genomic stability research.

[Learn More](#)

| Parameter        | Value                                                                                                           |
|------------------|-----------------------------------------------------------------------------------------------------------------|
| Product Name     | [KO Validated] MSH2 Rabbit pAb                                                                                  |
| Remarks/Alias    | FCC1; COCA1; HNPCC; LCF52; MSH-2; hMSH2; HNPCC1; LYNCH1; MMRC52; H2                                             |
| Species          | Human                                                                                                           |
| GeneID (human)   | 4436                                                                                                            |
| GeneID           | 4436                                                                                                            |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of human MSH2 (NP_000242.1) |
| Source           | Rabbit                                                                                                          |
| Category         | Polyclonal Antibodies                                                                                           |
| Application      | WB, IF/ICC, IP, ELISA                                                                                           |
| Cross Reactivity | Human, Rat                                                                                                      |
| SWISS            | P43246                                                                                                          |
| Protein Weight   | 105kDa                                                                                                          |
| Shipping         | Ice bag                                                                                                         |