

Anti-Msh4 Rabbit Polyclonal Antibody For Wb And Elisa - O15457

Item Number: CM0025133

Introduction

Rabbit polyclonal MSH4 antibody for WB and ELISA, reactive with human, mouse, and rat. MSH4 is essential for meiotic recombination and chromosome segregation, making it a key research tool for meiosis and fertility studies.

[Learn More](#)

Parameter	Value
Product Name	MSH4 Rabbit pAb
Remarks/Alias	ASG; POF20; SPGF2; MSH4
Species	Human
Gene ID	4438
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 697-936 of human MSH4 (NP_002431.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	O15457
Protein Weight	105kDa
Shipping	Ice bag