

## Anti-Msh6 Monoclonal Antibody For Wb, If/Icc, Elisa - P52701

Item Number: CM0011429



### Introduction

Anti-MSH6 Rabbit monoclonal antibody, validated for WB, IF/ICC, ELISA in human, mouse, rat; targets DNA mismatch repair protein Msh6 (P52701), ideal for cancer research and diagnostics.

[Learn More](#)

| Parameter        | Value                                                                                                                    |
|------------------|--------------------------------------------------------------------------------------------------------------------------|
| Product Name     | MSH6 Rabbit mAb                                                                                                          |
| Remarks/Alias    | GTBP; HSAP; p160; GTMBP; MSH-6; HNPCC5; LYNCH5; MMRC53; MSH6                                                             |
| Species          | Human                                                                                                                    |
| GeneID           | 2956                                                                                                                     |
| Immunogen        | Synthetic Peptide - A synthetic peptide corresponding to a sequence within amino acids 1-100 of human MSH6 (NP_000170.1) |
| Source           | Rabbit                                                                                                                   |
| Category         | Monoclonal Antibodies                                                                                                    |
| Application      | WB, IF/ICC, ELISA                                                                                                        |
| Cross Reactivity | Human, Mouse, Rat                                                                                                        |
| SWISS            | P52701                                                                                                                   |
| Protein Weight   | 153kDa                                                                                                                   |
| Shipping         | Ice bag                                                                                                                  |