

Anti-Msh6 Rabbit Monoclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P52701

Item Number: CM0006364



Introduction

KO Validated MSH6 Rabbit Monoclonal Antibody for DNA mismatch repair research. Applications: WB, IHC-P, IF/ICC, ELISA. Reacts with Human, Mouse, Rat. ~153kDa.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] MSH6 Rabbit mAb
Remarks/Alias	GTBP; HSAP; p160; GTMBP; MSH-6; HNPCC5; LYNCH5; MMRC53
Species	Human
GeneID (Human)	2956
GeneID	2956
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human MSH6 (NP_000170.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P52701
Protein Weight	153kDa
Shipping	Ice bag