

Anti-Msh6 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P52701

Item Number: CM0021938



Introduction

Rabbit anti-MSH6 polyclonal antibody validated for Western blot, IHC-P, IF/ICC, and ELISA. Detects human, mouse, and rat MSH6 (DNA mismatch repair protein). Suitable for MMR research.

[Learn More](#)

Parameter	Value
Product Name	MSH6 Rabbit pAb
Remarks/Alias	GTBP; HSAP; p160; GTMBP; MSH-6; HNPCC5; LYNCH5; MMRC53; MSH6
Species	Human
Gene ID (Human)	2956
Gene ID	2956
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-290 of human MSH6 (NP_000170.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P52701
Protein Weight	153kDa
Shipping	Ice bag