

Anti-Msh6 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P52701

Item Number: CM0020670



Introduction

High-quality rabbit polyclonal antibody against MSH6 for WB, IHC-P, IF/ICC, and ELISA applications. Detects MSH6 in human, mouse, rat, and monkey samples. A reliable tool for DNA mismatch repair research.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] MSH6 Rabbit pAb
Remarks/Alias	GTBP; HSAP; p160; GTMBP; MSH-6; HNPCC5; LYNCH5; MMRCS3; MSH6
Species	Human
GeneID (Human)	2956
GeneID	2956
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-100 of human MSH6 (NP_000170.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat, MK
SWISS	P52701
Protein Weight	153kDa
Shipping	Ice bag