

Anti-Msin3A Monoclonal Antibody For Ihc-P, If/Icc, Elisa - Q96St3

Item Number: CM0007315

Introduction

High-quality rabbit monoclonal antibody targeting human mSin3A (SWISS Q96ST3). Suitable for IHC-P, IF/ICC, and ELISA. Cross-reacts with mouse and rat. Ideal for transcriptional repression and circadian rhythm research.

[Learn More](#)

Parameter	Value
Product Name	mSin3A Rabbit mAb
Remarks/Alias	WITKOS; DEL15Q24; CHR15DELq24; mSin3A
Species	Human
Gene ID (Human)	25942
Gene ID	25942
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 1174-1273 of human mSin3A (Q96ST3).
Source	Rabbit
Category	Monoclonal Antibodies
Application	IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q96ST3
Protein Weight	145kDa
Shipping	Ice bag