

Anti-Msl1 Rabbit Polyclonal Antibody For Wb, Elisa - Q68Dk7

Item Number: CM0020973



Introduction

High-quality rabbit polyclonal antibody against human MSL1, validated for Western blot and ELISA. Cross-reacts with mouse and rat. Ideal for studying histone acetylation, chromatin stability, and gene dosage compensation pathways.

[Learn More](#)

Parameter	Value
Product Name	MSL1 Rabbit pAb
Remarks/Alias	MSL-1; MSL1
Species	Human
Gene ID (Human)	339287
Gene ID	339287
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 210-300 of human MSL1 (NP_001012241)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q68DK7
Protein Weight	67kDa
Shipping	Ice bag