

Anti-Msmb Rabbit Monoclonal Antibody For Wb, If/Icc, Elisa - P08118

Item Number: CM0008600



Introduction

High-quality rabbit monoclonal antibody against human beta-microseminoprotein (MSMB), validated for WB, IF/ICC, and ELISA. Cross-reacts with mouse MSMB. Ideal for prostate and seminal plasma research.

[Learn More](#)

Parameter	Specification
Product Name	MSMB Rabbit mAb
Remarks/Alias	MSP; PSP; IGBF; MSPB; PN44; PRPS; HPC13; PSP57; PSP94; PSP-94; MSMB
Species	Human
Gene ID	4477
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 1-100 of human MSMB (P08118).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	P08118
Protein Weight	13kDa
Shipping	Ice bag