


Anti-Msp/Mst1 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P26927

Item Number: CM0016863



Introduction

Rabbit polyclonal antibody against human MSP/MST1 (Hepatocyte growth factor-like protein), validated for WB, IHC-P, IF/ICC, ELISA. Cross-reacts with mouse and rat. Ideal for studying macrophage stimulation and liver regeneration.

[Learn More](#)

Parameter	Value
Product Name	MSP/MST1 Rabbit pAb
Remarks/Alias	MSP; HGFL; NF15S2; D3F15S2; DNF15S2; MSP/MST1
Species	Human
GeneID (Human)	4485
GeneID	4485
Immunogen	Recombinant fusion protein containing amino acids 200-500 of human MSP/MST1 (NP_066278.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P26927
Protein Weight	80kDa
Shipping	Ice bag