


## Anti-Msr1 Rabbit Polyclonal Antibody For Wb And Elisa - P21757

Item Number: CM0017796



### Introduction

Rabbit polyclonal antibody against human MSR1, suitable for WB and ELISA. Cross-reacts with mouse and rat. Detects a 50 kDa protein. Useful for studying macrophage scavenger receptor function in atherosclerosis and innate immunity.

[Learn More](#)

| Parameter        | Value                                                                                                             |
|------------------|-------------------------------------------------------------------------------------------------------------------|
| Product Name     | MSR1 Rabbit pAb                                                                                                   |
| Remarks/Alias    | SRA; SR-A; CD204; SR-AI; phSR1; phSR2; SCARA1; SR-AII; SR-AIII; MSR1                                              |
| Species          | Human                                                                                                             |
| GeneID           | 4481                                                                                                              |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 192-451 of human MSR1 (NP_619729.1) |
| Source           | Rabbit                                                                                                            |
| Category         | Polyclonal Antibodies                                                                                             |
| Application      | WB, ELISA                                                                                                         |
| Cross Reactivity | Human, Mouse, Rat                                                                                                 |
| SWISS            | P21757                                                                                                            |
| Protein Weight   | 50 kDa                                                                                                            |
| Shipping         | Ice bag                                                                                                           |