

Anti-Msra (Pmsr) Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q9Uj68

Item Number: CM0014420

Introduction

Rabbit polyclonal antibody against MSRA (methionine sulfoxide reductase A). Validated for WB, IF/ICC, ELISA. Cross-reacts with human, mouse, rat. Targets a mitochondrial and cytosolic repair enzyme critical for oxidative stress protection. 26 kDa.

[Learn More](#)

Parameter	Value
Product Name	MSRA Rabbit pAb
Remarks/Alias	PMSR; MSRA
Species	Human
GeneID (Human)	4482
GeneID	4482
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 24-235 of human MSRA (NP_036463.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UJ68
Protein Weight	26kDa
Shipping	Ice bag