

Anti-Msra Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q9Uj68

Item Number: CM0025386

Introduction

High-quality rabbit polyclonal antibody targeting human MSRA (PMSR), validated for WB, IHC-P, IF/ICC, and ELISA. Cross-reacts with human, mouse, rat. Ideal for oxidative stress and protein repair research.

[Learn More](#)

Parameter	Value
Product Name	MSRA Rabbit pAb
Remarks/Alias	PMSR; MSRA
Species	Human
Gene ID (Human)	4482
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 24-235 of human MSRA (NP_036463.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UJ68
Protein Weight	26kDa
Shipping	Ice bag