

Anti-Msrb1 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q9Nzv6

Item Number: CM0027858



Introduction

Rabbit polyclonal antibody against human MSRB1 (Selenoprotein X). Validated for WB, IHC-P, IF/ICC, ELISA applications. Reacts with human, mouse, and rat MSRB1. Suitable for research in oxidative stress, actin dynamics, and innate immunity.

[Learn More](#)

Parameter	Value
Product Name	MSRB1 Rabbit pAb
Remarks/Alias	SELR; SELX; SepR; SEPX1; HSPC270; SELENOR; SELENOX; MSRB1
Species	Human
GeneID (Human)	51734
GeneID	51734
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-94 of human MSRB1 (NP_057416.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9NZV6
Protein Weight	13kDa
Shipping	Ice bag