

Anti-Msrb2 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q9Y3D2

Item Number: CM0026426

Introduction

Rabbit polyclonal antibody against human MSRB2, a mitochondrial methionine-R-sulfoxide reductase involved in oxidative stress protection. Validated for WB, IHC-P, ELISA; cross-reacts with mouse and rat.

[Learn More](#)

Parameter	Value
Product Name	MSRB2 Rabbit pAb
Remarks/Alias	CBS1; MSRB; PILB; CBS-1; CGI-131; MSRB2
Species	Human
Gene ID (Human)	22921
Gene ID	22921
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 10-182 of human MSRB2 (NP_036360.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y3D2
Protein Weight	20kDa
Shipping	Ice bag