


Anti-MsrB3 Rabbit Monoclonal Antibody For Wb And Elisa - Q8IXL7

Item Number: CM0002898



Introduction

MSRB3 Rabbit monoclonal antibody for WB and ELISA. Cross-reacts with human, mouse, rat. Ideal for studies on oxidative stress, methionine reduction, and hearing loss.

[Learn More](#)

Parameter	Value
Product Name	MSRB3 Rabbit mAb
Remarks/Alias	DFNB74; MSRB3
Species	Human
GeneID	253827
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 50-192 of human MSRB3 (Q8IXL7)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8IXL7
Protein Weight	21kDa
Shipping	Ice bag