


Anti-Msrb3 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q8IXL7

Item Number: CM0026814



Introduction

Rabbit polyclonal antibody targeting MSRB3 (Methionine-R-sulfoxide reductase B3) for WB, IF/ICC, ELISA applications. Cross-reacts with human, mouse, and rat. Suitable for oxidative stress and deafness research.

[Learn More](#)

Parameter	Value
Product Name	MSRB3 Rabbit pAb
Remarks/Alias	DFNB74; MSRB3
Species	Human
GeneID (Human)	253827
GeneID	253827
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-185 of human MSRB3 (NP_001026849.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8IXL7
Protein Weight	21kDa
Shipping	Ice bag