


Anti-Mst1 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q13043

Item Number: CM0023708



Introduction

Rabbit polyclonal anti-MST1 (STK4) antibody, validated for WB, IHC-P, IF/ICC, ELISA. Reacts with human, mouse, rat. Ideal for studying Hippo signaling, apoptosis, and cancer research.

[Learn More](#)

Parameter	Value
Product Name	MST1 Rabbit pAb
Remarks/Alias	KRS2; MST1; YSK3; STK4
Species	Human
Gene ID (Human)	6789
Gene ID	6789
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 303-487 of human MST1 (NP_006273.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13043
Protein Weight	56kDa
Shipping	Ice bag