


Anti-Mst1/Stk4 Polyclonal Antibody For Wb, Elisa - Q13043

Item Number: CM0023659



Introduction

KO Validated rabbit polyclonal antibody recognizing human MST1/STK4 (KRS2, MST1, YSK3). Confirmed for WB and ELISA; cross-reacts with human, mouse, and rat. For Hippo signaling and apoptosis studies.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] MST1/STK4 Rabbit pAb
Remarks/Alias	KRS2; MST1; YSK3; K4
Species	Human
Gene ID (Human)	6789
Gene ID	6789
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 303-487 of human MST1/STK4 (NP_006273.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13043
Protein Weight	56kDa
Shipping	Ice bag