


## Anti-Mst1/Stk4 Rabbit Polyclonal Antibody For Wb, If/Icc, Ip, Elisa - Q13043

Item Number: CM0029899



### Introduction

High-quality rabbit polyclonal antibody against MST1/STK4, validated for WB, IF/ICC, IP, ELISA in human, mouse, rat. Ideal for apoptosis and Hippo signaling research.

[Learn More](#)

| Parameter        | Value                                                                                                                   |
|------------------|-------------------------------------------------------------------------------------------------------------------------|
| Product Name     | MST1/STK4 Rabbit pAb                                                                                                    |
| Remarks/Alias    | KRS2; MST1; YSK3; K4                                                                                                    |
| Species          | Human                                                                                                                   |
| GeneID (Human)   | 6789                                                                                                                    |
| GeneID           | 6789                                                                                                                    |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 303-487 of human MST1/STK4 (NP_006273.1). |
| Source           | Rabbit                                                                                                                  |
| Category         | Polyclonal Antibodies                                                                                                   |
| Application      | WB, IF/ICC, IP, ELISA                                                                                                   |
| Cross Reactivity | Human, Mouse, Rat                                                                                                       |
| SWISS            | Q13043                                                                                                                  |
| Protein Weight   | 56 kDa                                                                                                                  |
| Shipping         | Ice bag                                                                                                                 |