


Anti-Msy2 Rabbit Polyclonal Antibody For Wb And Elisa - Q9Y2T7

Item Number: CM0019814



Introduction

Rabbit polyclonal antibody targeting human MSY2 (Y-box-binding protein 2, DBPC, CSDA3, CONTRIN). Validated for Western blot and ELISA; cross-reacts with human and mouse. Detects endogenous MSY2 at ~39 kDa. Ideal for germ cell biology and cancer research.

[Learn More](#)

Parameter	Value
Product Name	MSY2 Rabbit pAb
Remarks/Alias	DBPC; MSY2; CSDA3; CONTRIN
Species	Human
Gene ID	51087
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-100 of human MSY2 (NP_057066.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9Y2T7
Protein Weight	39kDa
Shipping	Ice bag