

Anti-Mt-Atp6 Polyclonal Antibody For Wb, If/Icc, Ihc-P, Elisa - P00846

Item Number: CM0021705



Introduction

High-quality rabbit polyclonal antibody against human MT-ATP6 (ATP synthase F(0) subunit a). Applicable in WB, IF/ICC, IHC-P, ELISA. Cross-reacts with mouse and rat. Suitable for mitochondrial disease research.

[Learn More](#)

Parameter	Value
Product Name	MT-ATP6 Rabbit pAb
Remarks/Alias	ATPase6; MTATP6; ATP6; MT-ATP6
Species	Human
Gene ID (Human)	4508
Gene ID	4508
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human MT-ATP6 (YP_003024031.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P00846
Protein Weight	25 kDa
Shipping	Ice bag