


Anti-Mt-Co3 Rabbit Pab For Wb, Ihc-P, Elisa - P00414

Item Number: CM0022567



Introduction

Rabbit polyclonal antibody targeting MT-CO3 (Cytochrome c oxidase subunit 3), suitable for WB, IHC-P, ELISA in human, mouse, and rat samples. Ideal for mitochondrial respiratory chain research.

[Learn More](#)

Parameter	Value
Product Name	MT-CO3 Rabbit pAb
Remarks/Alias	COIII; MTCO3; COX3; MT-CO3
Species	Human
Gene ID (Human)	4514
Gene ID	4514
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human MT-CO3 (YP_003024032.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P00414
Protein Weight	30kDa
Shipping	Ice bag