


Anti-Mt-Cyb Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P00156

Item Number: CM0024123



Introduction

Rabbit polyclonal antibody against human Cytochrome b (MT-CYB), a mitochondrial respiratory chain complex III component. Suitable for WB, IHC-P, IF/ICC, ELISA. Reacts with human, mouse, rat. Detects ~43 kDa band.

[Learn More](#)

Parameter	Value
Product Name	MT-CYB Rabbit pAb
Remarks/Alias	MTCYB; CYTB; MT-CYB
Species	Human
Gene ID	4519
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 100-200 of human MT-CYB (YP_003024038.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P00156
Protein Weight	43kDa
Shipping	Ice bag