


Anti-Mt-Cyb Rabbit Polyclonal Antibody For Wb And Elisa - P00158

Item Number: CM0025409



Introduction

Rabbit polyclonal antibody against mouse MT-CYB (Cytochrome b), validated for WB and ELISA, cross-reacts with human and rat. Ideal for mitochondrial respiratory chain and oxidative phosphorylation research.

[Learn More](#)

Parameter	Value
Product Name	MT-CYB Rabbit pAb
Remarks/Alias	CYTB; MT-CYB
Species	Mouse
GeneID (Human)	4519
GeneID	17711
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 3-103 of mouse MT-CYB (NP_904340.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P00158
Protein Weight	43kDa
Shipping	Ice bag