


## Anti-Mt-Nd1 Rabbit Monoclonal Antibody (Wb, Elisa) - P03886

Item Number: CM0005204



### Introduction

Rabbit monoclonal antibody against human MT-ND1, validated for Western blot and ELISA. Cross-reacts with mouse and rat. Suitable for studying mitochondrial complex I function, protein weight ~36kDa.

[Learn More](#)

| Parameter        | Value                                                                                             |
|------------------|---------------------------------------------------------------------------------------------------|
| Product Name     | MT-ND1 Rabbit mAb                                                                                 |
| Remarks/Alias    | MTND1, ND1, MT-ND1                                                                                |
| Species          | Human                                                                                             |
| GeneID (human)   | 4535                                                                                              |
| GeneID           | 4535                                                                                              |
| Immunogen        | A synthetic peptide corresponding to a sequence within amino acids 1-100 of human MT-ND1 (P03886) |
| Source           | Rabbit                                                                                            |
| Category         | Monoclonal Antibodies                                                                             |
| Application      | WB, ELISA                                                                                         |
| Cross Reactivity | Mouse, Rat                                                                                        |
| SWISS            | P03886                                                                                            |
| Protein Weight   | 36kDa                                                                                             |
| Shipping         | Ice bag                                                                                           |