


Anti-Mt-Nd2 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P03891

Item Number: CM0023453



Introduction

Rabbit polyclonal antibody targeting human MT-ND2 (NADH dehydrogenase subunit 2), a core Complex I subunit. Validated for WB, IHC-P, IF/ICC, ELISA. Cross-reacts with mouse and rat.

[Learn More](#)

Parameter	Value
Product Name	MT-ND2 Rabbit pAb
Remarks/Alias	MTND2; ND2; MT-ND2
Species	Human
Gene ID (Human)	4536
Gene ID	4536
Immunogen	Synthetic peptide
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P03891
Protein Weight	39kDa
Shipping	Ice bag