


Anti-Mt-Nd3 Rabbit Polyclonal Antibody For Wb - P03897

Item Number: CM0023992



Introduction

Rabbit polyclonal antibody targeting human MT-ND3, a core mitochondrial Complex I subunit. Validated for WB, IHC-P, IF/ICC, IP, ELISA. Cross-reacts with mouse and rat. Suitable for studying mitochondrial function and respiratory chain disorders.

[Learn More](#)

Parameter	Value
Product Name	MT-ND3 Rabbit pAb
Remarks/Alias	MTND3; ND3; MT-ND3
Species	Human
GeneID (Human)	4537
GeneID	4537
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-70 of human MT-ND3 (YP_003024033.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P03897
Protein Weight	13kDa
Shipping	Ice bag